

## Prezalumab

## Chemical Properties

CAS No. :	1523164-68-2
Formula:	
Molecular Weight:	145.42 kDa
Storage:	Store at low temperature -20°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Prezalumab (AMG 557) is a human IgG2 monoclonal antibody targeting ICOSL and BAFF, used for studying Sjogren's syndrome, cutaneous lupus erythematosus, psoriasis, systemic lupus erythematosus arthritis, and systemic lupus erythematosus (SLE).
Targets(IC50)	Others, Interleukin
In vivo	In mouse NZB/NZW lupus model and collagen-induced arthritis (CIA) models, Prezalumab demonstrates more efficacious effects than a single ICOSL or BAFF inhibitor [1]. Prezalumab, similar to the pharmacological effect of BAFF inhibition, reduces B cells in cynomolgus monkeys[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0069 mL	0.0344 mL	0.0688 mL
5 mM	0.0014 mL	0.0069 mL	0.0138 mL
10 mM	0.0007 mL	0.0034 mL	0.0069 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481