


## Siraitic Acid B

## Chemical Properties

CAS No. :	183374-16-5
Formula:	C <sub>29</sub> H <sub>42</sub> O <sub>5</sub>
Molecular Weight:	470.64
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Siraitic Acid B is a natural product(cucurbitane triterpenoid) isolated from the root of <i>S. grosvenori</i> .
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1248 mL	10.6238 mL	21.2477 mL
5 mM	0.425 mL	2.1248 mL	4.2495 mL
10 mM	0.2125 mL	1.0624 mL	2.1248 mL
50 mM	0.0425 mL	0.2125 mL	0.425 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481