


Tyrosol

Chemical Properties

CAS No. :	501-94-0
Formula:	C ₈ H ₁₀ O ₂
Molecular Weight:	138.16
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tyrosol (2-(4-Hydroxyphenyl)ethanol) is a phenolic antioxidant that is naturally present in several foods such as wines and green tea and is present most abundantly in olives.
Targets(IC50)	Antioxidant, NF-κB, Endogenous Metabolite

Solubility Information

Solubility	DMSO: 255 mg/mL (1845.69 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (14.48 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.238 mL	36.1899 mL	72.3798 mL
5 mM	1.4476 mL	7.238 mL	14.476 mL
10 mM	0.7238 mL	3.619 mL	7.238 mL
50 mM	0.1448 mL	0.7238 mL	1.4476 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Chandramohan R, Pari L, Rathinam A, Sheikh BA. Chem Biol Interact. 2015 Mar 5;229:44-54.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481