

(2E)-4-[4,5-bis(hydroxymethyl)-3-methoxy-2-methylphenoxy]-2-methylbut-2-en-1-ol

Chemical Properties


CAS No. :

Formula: C₁₅H₂₂O₅

Molecular Weight: 282.34

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

(2E)-4-[4,5-bis(hydroxymethyl)-3-methoxy-2-methylphenoxy]-2-methylbut-2-en-1-ol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5418 mL	17.7091 mL	35.4183 mL
5 mM	0.7084 mL	3.5418 mL	7.0837 mL
10 mM	0.3542 mL	1.7709 mL	3.5418 mL
50 mM	0.0708 mL	0.3542 mL	0.7084 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481