

Compound T128663

Chemical Properties

CAS No. :

Formula: C₃₅H₅₈O₁₀

Molecular Weight: 638.84

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	3-[3-(5-hydroperoxy-6-methyl-2-(((2S,3R,4S,5S)-3,4,5-trihydroxyoxan-2-yl)oxy)hept-6-en-2-yl)-4-hydroxy-6,9a,9b-trimethyl-7-(prop-1-en-2-yl)-dodecahydro-1H-cyclopenta[a]naphthalen-6-yl]propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5653 mL	7.8267 mL	15.6534 mL
5 mM	0.3131 mL	1.5653 mL	3.1307 mL
10 mM	0.1565 mL	0.7827 mL	1.5653 mL
50 mM	0.0313 mL	0.1565 mL	0.3131 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481