


## Butyrophenone

## Chemical Properties

CAS No. :	495-40-9
Formula:	C <sub>10</sub> H <sub>12</sub> O
Molecular Weight:	148.2
Storage:	Pure form: -20°C for 3 years   In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Butyrophenone (NSC-8463) is a chemical compound. Some of its derivatives are used to treat various psychiatric disorders.
Targets(IC50)	Others,Dopamine Receptor

## Solubility Information

Solubility	DMSO: 50 mg/mL (337.38 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.7476 mL	33.7382 mL	67.4764 mL
5 mM	1.3495 mL	6.7476 mL	13.4953 mL
10 mM	0.6748 mL	3.3738 mL	6.7476 mL
50 mM	0.135 mL	0.6748 mL	1.3495 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481