


Simiarenone

Chemical Properties

CAS No. :	2318-78-7
Formula:	C ₃₀ H ₄₈ O
Molecular Weight:	424.7
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Simiarenone is a natural product of Euphorbia, Euphorbiaceae. The catalog number is TN5025 and the CAS number is 2318-78-7. Simiarenone can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3546 mL	11.773 mL	23.546 mL
5 mM	0.4709 mL	2.3546 mL	4.7092 mL
10 mM	0.2355 mL	1.1773 mL	2.3546 mL
50 mM	0.0471 mL	0.2355 mL	0.4709 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481