


5-(Prop-1-ynyl)pyridine-3-ylboronic acid

Chemical Properties

CAS No. :	917471-30-8
Formula:	C ₈ H ₈ BNO ₂
Molecular Weight:	160.97
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Intermediate and building block-boronic acid, heterocyclic compound-pyridine
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.2123 mL	31.0617 mL	62.1234 mL
5 mM	1.2425 mL	6.2123 mL	12.4247 mL
10 mM	0.6212 mL	3.1062 mL	6.2123 mL
50 mM	0.1242 mL	0.6212 mL	1.2425 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481