

Sanggenon O

Chemical Properties


CAS No. :

Formula: C₄₀H₃₆O₁₂

Molecular Weight: 708.716

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sanggenon O is a useful organic compound for research related to life sciences and the catalog number is T126526.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.411 mL	7.055 mL	14.1099 mL
5 mM	0.2822 mL	1.411 mL	2.822 mL
10 mM	0.1411 mL	0.7055 mL	1.411 mL
50 mM	0.0282 mL	0.1411 mL	0.2822 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481