

## Sanggenon O

## Chemical Properties

CAS No. :	101664-32-8
Formula:	C40H36O12
Molecular Weight:	708.71
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Sanggenon O is a useful organic compound for research related to life sciences and the catalog number is T126526.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.411 mL	7.0551 mL	14.1101 mL
5 mM	0.2822 mL	1.411 mL	2.822 mL
10 mM	0.1411 mL	0.7055 mL	1.411 mL
50 mM	0.0282 mL	0.1411 mL	0.2822 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481