


Loprodiol

Chemical Properties

CAS No. : 2209-86-1
 Formula: C₅H₁₀Cl₂O₂
 Molecular Weight: 173.04
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Loprodiol is a biochemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.779 mL	28.8951 mL	57.7901 mL
5 mM	1.1558 mL	5.779 mL	11.558 mL
10 mM	0.5779 mL	2.8895 mL	5.779 mL
50 mM	0.1156 mL	0.5779 mL	1.1558 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481