

Compound T130573

Chemical Properties


CAS No. :

Formula: C₂₆H₄₂O₈

Molecular Weight: 482.61

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,4a,7-trimethyl-7-(2-([3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)ethyl)-1,2,3,4,4a,4b,5,6,7,8,10,10a-dodecahydrophenanthrene-1-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0721 mL	10.3603 mL	20.7207 mL
5 mM	0.4144 mL	2.0721 mL	4.1441 mL
10 mM	0.2072 mL	1.036 mL	2.0721 mL
50 mM	0.0414 mL	0.2072 mL	0.4144 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481