


GW300657X

Chemical Properties

CAS No. : 388626-82-2
 Formula: C₂₁H₁₈N₆O₄S
 Molecular Weight: 450.47
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	GW300657X is a small molecule inhibitor that blocks toxoplasma gondii rhoptyr kinase 18 and maintains moderate CDK2-cyclin A inhibitory activity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2199 mL	11.0995 mL	22.199 mL
5 mM	0.444 mL	2.2199 mL	4.4398 mL
10 mM	0.222 mL	1.110 mL	2.2199 mL
50 mM	0.0444 mL	0.222 mL	0.444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481