


## 5-Phenylcytidine

## Chemical Properties

CAS No. :	83866-19-7
Formula:	C <sub>15</sub> H <sub>17</sub> N <sub>3</sub> O <sub>5</sub>
Molecular Weight:	319.31
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	5-Phenylcytidine is a Nucleoside Derivative - 5-Modified pyrimidine nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1318 mL	15.6588 mL	31.3175 mL
5 mM	0.6264 mL	3.1318 mL	6.2635 mL
10 mM	0.3132 mL	1.5659 mL	3.1318 mL
50 mM	0.0626 mL	0.3132 mL	0.6264 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481