

1A-116

## Chemical Properties


CAS No. : 1430208-73-3

Formula: C<sub>16</sub>H<sub>16</sub>F<sub>3</sub>N<sub>3</sub>

Molecular Weight: 307.31

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1A-116 is a specific Rac1 inhibitor.
Targets(IC50)	Apoptosis,Rho,Ras
In vitro	1A-116 shows lesser effect on MCF7::pcDNA.3 cells than on MCF7::C1199 cells. 1A-116 treatment decreases phospho-PAK1 levels in a time-dependent manner. The presence of 1A-116 reverts the PAK1 phosphorylation induced by 4-hydroxytamoxifen (Tam). The presence of 1A-116 also effectively reverts Rac1-PAK1-mediated estrogen receptor (ER) phosphorylation at Ser305[1]. 1A-116 shows a significant increase in antiproliferative activity on F3II cells, showing an IC50 value of 4 μM. A-116 also dramatically impairs Rac1 activation at low micromolar range (1 μM)[2].
In vivo	Daily treatment of mice with compound 1A-116 at 3mg/kg body weight/day reduces about 60% the formation of total metastatic lung colonies. A significant antitumor activity is obtained for macronodules (more than 1 mm in diameter) by treatment with 1A-116 in this highly aggressive breast cancer model. The treatment with 1A-116 reduces the total lung weight compare to the control group, leading to a total weight similar to the average pulmonary weight of Balb/c mice[2].

## Solubility Information

Solubility	DMSO: 100 mg/mL (325.4 mM),Sonication is recommended. H2O: < 0.1 mg/mL (insoluble) (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 4 mg/mL (13.02 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

### Preparing Stock Solutions

---

	1mg	5mg	10mg
1 mM	3.254 mL	16.2702 mL	32.5404 mL
5 mM	0.6508 mL	3.254 mL	6.5081 mL
10 mM	0.3254 mL	1.627 mL	3.254 mL
50 mM	0.0651 mL	0.3254 mL	0.6508 mL

---

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

### Reference

Gonzalez N, et al. Pharmacological inhibition of Rac1-PAK1 axis restores tamoxifen sensitivity in human resistant breast cancer cells. *Cell Signal*. 2017 Jan;30:154-161.

Cardama GA, et al. Preclinical development of novel Rac1-GEF signaling inhibitors using a rational design approach in highly aggressive breast cancer cell lines. *Anticancer Agents Med Chem*. 2014;14(6):840-51.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481