


ICI-170809

Chemical Properties

CAS No. : 85275-48-5
 Formula: C₂₁H₂₄N₂S
 Molecular Weight: 336.49
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	ICI-170809 is a 5-HT ₂ receptor antagonist.
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9719 mL	14.8593 mL	29.7186 mL
5 mM	0.5944 mL	2.9719 mL	5.9437 mL
10 mM	0.2972 mL	1.4859 mL	2.9719 mL
50 mM	0.0594 mL	0.2972 mL	0.5944 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481