


IDX-989

Chemical Properties

CAS No. : 881924-55-6
 Formula: C₂₀H₁₆ClFN₃O₃P
 Molecular Weight: 431.78
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	IDX-989 is an HIV-1 non-nucleoside reverse transcriptase inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.316 mL	11.580 mL	23.1599 mL
5 mM	0.4632 mL	2.316 mL	4.632 mL
10 mM	0.2316 mL	1.158 mL	2.316 mL
50 mM	0.0463 mL	0.2316 mL	0.4632 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481