


WKYMVM

Chemical Properties

CAS No. :	187986-17-0
Formula:	C41H61N9O7S2
Molecular Weight:	856.11
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	WKYMVM is a N-formyl peptide receptor (FPR1) agonist.
Targets(IC50)	Apoptosis,Others,PERK
In vitro	Neutrophils are activated by WKYMVM through FPRL1 and that FPRL2 is a chemotactic receptor transducing signals in myeloid cells[1].
In vivo	Using blood samples withdrawn from 28 AL patients and 31 healthy controls, found that, in the absence of exogenous WKYMVM, WKYMVM receptor and bactericidal activity did not differ between patients and controls.Addition of WKYMVM markedly increased the bactericidal activities of these cells in a dose-dependent manner.During induction chemotherapy, there were significant increases in bactericidal activity in the presence and absence of 1nM WKYMVM, with higher bactericidal activities at the time of complete remission than at the time of diagnosis or on day 15.During consolidation chemotherapy, WKYMVM had no effect on bactericidal activities.Patients showed significant increase in the concentrations of TNF alpha, IL-1b, IL-6 and IL-8, but significant decrease in the concentrations of IL-2, IL-4 and IL-12.TNF alpha, IL-1b and IL-6 showed significant negative correlations with bactericidal activities of patient neutrophils at time of diagnosis, and IL-4 showed a significant positive correlation with bactericidal activities,these findings indicate that WKYMVM enhances bactericidal activity in patients with AL[2].
Cell Research	Ligand-mediated internalization of FPRL1 and FPRL2 was followed by an indirect immunofluorescence microscopy staining technique.?Cells were seeded on polyornithine-coated coverslips 48 h prior to the experiments.?To stop protein synthesis, cells were incubated in RPMI 1640 medium containing 1% BSA and 100 µg/ml cycloheximide for 3 h and maintained in the same medium throughout the experiment.? To study receptor internalization, RINm5F-FPRL1 cells and RINm5F-FPRL2 cells were treated with either WKYMVM or WKYMVm at 37 °C for various time periods ranging from 0 to 30 min. Cells were finally fixed and permeabilized with acetone at ?20 °C for 30 s[1].

Solubility Information

A DRUG SCREENING EXPERT

Solubility	DMSO: 10 mM, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1681 mL	5.8404 mL	11.6807 mL
5 mM	0.2336 mL	1.1681 mL	2.3361 mL
10 mM	0.1168 mL	0.584 mL	1.1681 mL
50 mM	0.0234 mL	0.1168 mL	0.2336 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Christophe T , Karlsson A , Dugave C , et al. The Synthetic Peptide Trp-Lys-Tyr-Met-Val-Met-NH₂ Specifically Activates Neutrophils through FPRL1/Lipoxin A4 Receptors and Is an Agonist for the Orphan Monocyte-expressed Chemoattractant Receptor FPRL2[J]. Journal of Biological Chemistry, 2001, 276(24):21585-21593.

Kim H , Noh E K , Lee E J , et al. Enhanced bactericidal function by WKYMVm in patients with acute leukemia[J]. Leukemia Research, 2008, 32(5):717-725.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481