


RP-52770

Chemical Properties

CAS No. : 93363-02-1
 Formula: C₁₈H₁₄ClN₃O₂
 Molecular Weight: 355.84
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|-----------------------------------|
| Description | RP-52770 is a bioactive chemical. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.8103 mL | 14.0513 mL | 28.1025 mL |
| 5 mM | 0.5621 mL | 2.8103 mL | 5.6205 mL |
| 10 mM | 0.281 mL | 1.4051 mL | 2.8103 mL |
| 50 mM | 0.0562 mL | 0.281 mL | 0.5621 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481