


## Octenidine saccharin

## Chemical Properties

CAS No. : 86767-75-1  
 Formula: C<sub>50</sub>H<sub>72</sub>N<sub>6</sub>O<sub>6</sub>S<sub>2</sub>  
 Molecular Weight: 917.27  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Octenidine saccharin is a Dental Plaque Inhibitor
Targets(IC50)	Others,Antibacterial

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0902 mL	5.451 mL	10.9019 mL
5 mM	0.218 mL	1.0902 mL	2.1804 mL
10 mM	0.109 mL	0.5451 mL	1.0902 mL
50 mM	0.0218 mL	0.109 mL	0.218 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481