

PCC0208018

## Chemical Properties


CAS No. : 1673534-73-0

Formula: C7H13N5O3

Molecular Weight: 215.21

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PCC0208018 is a novel activator of effector T cells, enhancing T cell proliferation and activation to release interferon gamma (IFN- $\gamma$ ) and interleukin-2 (IL-2) without blocking the programmed cell death 1 (PD-1)/programmed cell death-ligand 1 (PD-L1) binding and not directly affecting tumor cell viability in vitro.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6466 mL	23.2331 mL	46.4662 mL
5 mM	0.9293 mL	4.6466 mL	9.2932 mL
10 mM	0.4647 mL	2.3233 mL	4.6466 mL
50 mM	0.0929 mL	0.4647 mL	0.9293 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481