

## Azido-C3-UV-biotin

## Chemical Properties

CAS No. :


Formula: C38H48N4O9S

Molecular Weight: 736.87

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Azido-C3-UV-biotin is a PEG-based linker for PROTACs, joining two essential ligands crucial for forming PROTAC molecules and enabling selective protein degradation via the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs comprise two distinct ligands linked together: one binds to an E3 ubiquitin ligase, and the other targets the specific protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3571 mL	6.7855 mL	13.5709 mL
5 mM	0.2714 mL	1.3571 mL	2.7142 mL
10 mM	0.1357 mL	0.6785 mL	1.3571 mL
50 mM	0.0271 mL	0.1357 mL	0.2714 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481