

Compound T130494

Chemical Properties


CAS No. :

Formula: C33H40O18

Molecular Weight: 724.67

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-(((2R,3S,4S,5R,6S)-3,4-dihydroxy-6-([5-hydroxy-2-(4-hydroxyphenyl)-4-oxo-3,4-dihydro-2H-1-benzopyran-7-yl]oxy)-5-(((2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy)oxan-2-yl]methoxy)-3-hydroxy-3-methyl-5-oxopentanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3799 mL	6.8997 mL	13.7994 mL
5 mM	0.276 mL	1.3799 mL	2.7599 mL
10 mM	0.138 mL	0.690 mL	1.3799 mL
50 mM	0.0276 mL	0.138 mL	0.276 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481