

Compound T127803

Chemical Properties

CAS No. :

Formula: C₂₂H₃₀O₁₁

Molecular Weight: 470.47

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(2S,3R,4S,5S,6R)-2-[[[(2R,3R,4S,5S,6R)-2-(deca-4,6,8-triyn-1-yloxy)-4,5-dihydroxy-6-(hydroxymethyl)oxan-3-yl]oxy]-6-(hydroxymethyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1255 mL	10.6277 mL	21.2553 mL
5 mM	0.4251 mL	2.1255 mL	4.2511 mL
10 mM	0.2126 mL	1.0628 mL	2.1255 mL
50 mM	0.0425 mL	0.2126 mL	0.4251 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481