

Pomalidomide-C4-OH

Chemical Properties

CAS No. :


Formula: C17H19N3O5

Molecular Weight: 345.36

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pomalidomide-C4-OH is a synthesized E3 ligase ligand-linker conjugate that incorporates the Pomalidomide based CRBN ligand and linker used in PROTAC technology.
Targets(IC50)	Others,E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8955 mL	14.4776 mL	28.9553 mL
5 mM	0.5791 mL	2.8955 mL	5.7911 mL
10 mM	0.2896 mL	1.4478 mL	2.8955 mL
50 mM	0.0579 mL	0.2896 mL	0.5791 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481