


## 2,4,6,6-Tetramethyl-3(6H)-pyridinone

## Chemical Properties

CAS No. :	203524-64-5
Formula:	C <sub>9</sub> H <sub>13</sub> NO
Molecular Weight:	151.21
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2,4,6,6-Tetramethyl-3(6H)-pyridinone, catalog number TN2709 and CAS number 203524-64-5, is a natural product used in life sciences research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.6133 mL	33.0666 mL	66.1332 mL
5 mM	1.3227 mL	6.6133 mL	13.2266 mL
10 mM	0.6613 mL	3.3067 mL	6.6133 mL
50 mM	0.1323 mL	0.6613 mL	1.3227 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481