


## 3(2H)-pyridazinone

## Chemical Properties

CAS No. :	504-30-3
Formula:	C <sub>4</sub> H <sub>4</sub> N <sub>2</sub> O
Molecular Weight:	96.089
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	3(2H)-pyridazinone, with CAS No. 504-30-3, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 3(2H)-pyridazinone provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	10.4069 mL	52.0346 mL	104.0691 mL
5 mM	2.0814 mL	10.4069 mL	20.8138 mL
10 mM	1.0407 mL	5.2035 mL	10.4069 mL
50 mM	0.2081 mL	1.0407 mL	2.0814 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481