

Compound 0080-0003

Chemical Properties


CAS No. :

Formula: C₁₉H₂₀O₅

Molecular Weight: 328.364

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound 0080-0003 is a useful organic compound for research related to life sciences and the catalog number is T131709.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0454 mL	15.2272 mL	30.4544 mL
5 mM	0.6091 mL	3.0454 mL	6.0909 mL
10 mM	0.3045 mL	1.5227 mL	3.0454 mL
50 mM	0.0609 mL	0.3045 mL	0.6091 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481