

## 3-Hydroxy-5-phenylpentanoic acid, O- $\beta$ -D-Glucopyran

### Chemical Properties


CAS No. :

Formula: C<sub>17</sub>H<sub>24</sub>O<sub>8</sub>

Molecular Weight: 356.371

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	3-Hydroxy-5-phenylpentanoic acid, O- $\beta$ -D-Glucopyran is a useful organic compound for research related to life sciences and the catalog number is T130838.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8061 mL	14.0304 mL	28.0607 mL
5 mM	0.5612 mL	2.8061 mL	5.6121 mL
10 mM	0.2806 mL	1.403 mL	2.8061 mL
50 mM	0.0561 mL	0.2806 mL	0.5612 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481