


Org 20599

## Chemical Properties

CAS No. : 156685-94-8  
 Formula: C<sub>25</sub>H<sub>40</sub>ClNO<sub>3</sub>  
 Molecular Weight: 438.04  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	GABAA receptor agonist
Targets(IC50)	Others,GABA Receptor

## Solubility Information

Solubility	DMSO: <43.8 mg/mL, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2829 mL	11.4145 mL	22.829 mL
5 mM	0.4566 mL	2.2829 mL	4.5658 mL
10 mM	0.2283 mL	1.1414 mL	2.2829 mL
50 mM	0.0457 mL	0.2283 mL	0.4566 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481