


Mal-NH-Boc

Chemical Properties

CAS No. :	134272-63-2
Formula:	C ₁₁ H ₁₆ N ₂ O ₄
Molecular Weight:	240.26
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Mal-NH-Boc is a alkyl/ether-based linker for PROTACs which joins two essential ligands, crucial for forming PROTAC molecules. This linker enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs are composed of two distinct ligands linked together: one targets an E3 ubiquitin ligase, and the other binds to a specific target protein. By leveraging the intracellular ubiquitin-proteasome system, PROTACs facilitate the selective degradation of target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1622 mL	20.8108 mL	41.6216 mL
5 mM	0.8324 mL	4.1622 mL	8.3243 mL
10 mM	0.4162 mL	2.0811 mL	4.1622 mL
50 mM	0.0832 mL	0.4162 mL	0.8324 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481