


## 3,5-Dimethyl-4-isoxazolecarbonitrile

## Chemical Properties

CAS No. :	31301-46-9
Formula:	C <sub>6</sub> H <sub>6</sub> N <sub>2</sub> O
Molecular Weight:	122.1246
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3,5-Dimethyl-4-isoxazolecarbonitrile, with CAS No. 31301-46-9, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 3,5-Dimethyl-4-isoxazolecarbonitrile provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.1887 mL	40.9433 mL	81.8867 mL
5 mM	1.6377 mL	8.1887 mL	16.3773 mL
10 mM	0.8189 mL	4.0943 mL	8.1887 mL
50 mM	0.1638 mL	0.8189 mL	1.6377 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481