

Compound N010-0010

Chemical Properties


CAS No. :

Formula: C₂₃H₃₈O₄

Molecular Weight: 378.553

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N010-0010 is a useful organic compound for research related to life sciences and the catalog number is T131365.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6417 mL	13.2083 mL	26.4166 mL
5 mM	0.5283 mL	2.6417 mL	5.2833 mL
10 mM	0.2642 mL	1.3208 mL	2.6417 mL
50 mM	0.0528 mL	0.2642 mL	0.5283 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481