


EI-1

## Chemical Properties

CAS No. : 301211-35-8  
 Formula: C<sub>16</sub>H<sub>19</sub>F<sub>3</sub>N<sub>6</sub>O<sub>3</sub>  
 Molecular Weight: 400.36  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	EI-1 is an inhibitor of Hepatitis C virus entry.
Targets(IC50)	Others,HCV Protease

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4978 mL	12.4888 mL	24.9775 mL
5 mM	0.4996 mL	2.4978 mL	4.9955 mL
10 mM	0.2498 mL	1.2489 mL	2.4978 mL
50 mM	0.050 mL	0.2498 mL	0.4996 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481