


Sulphenone

Chemical Properties

CAS No. :	80-00-2
Formula:	C ₁₂ H ₉ ClO ₂ S
Molecular Weight:	252.72
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sulphenone is an agent of Acaricide.
Targets(IC50)	Others,Parasite

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9569 mL	19.7847 mL	39.5695 mL
5 mM	0.7914 mL	3.9569 mL	7.9139 mL
10 mM	0.3957 mL	1.9785 mL	3.9569 mL
50 mM	0.0791 mL	0.3957 mL	0.7914 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481