

(2R,3R,4S,5S,6R)-2-[(8E)-dec-8-en-4,6-diyn-1-yloxy]-6-(hydroxymethyl)oxane-3,4,5-triol

Chemical Properties

CAS No. :

Formula: C₁₆H₂₂O₆

Molecular Weight: 310.35

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(2R,3R,4S,5S,6R)-2-[(8E)-dec-8-en-4,6-diyn-1-yloxy]-6-(hydroxymethyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2222 mL	16.1108 mL	32.2217 mL
5 mM	0.6444 mL	3.2222 mL	6.4443 mL
10 mM	0.3222 mL	1.6111 mL	3.2222 mL
50 mM	0.0644 mL	0.3222 mL	0.6444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481