

Compound T128142

Chemical Properties


CAS No. :

Formula: C₁₆H₂₂O₁₀S

Molecular Weight: 406.4

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

[(2S,3R,4S,5S,6R)-4,5-dihydroxy-6-(hydroxymethyl)-2-[2-methoxy-4-(prop-2-en-1-yl)phenoxy]oxan-3-yl]oxidanesulfonic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4606 mL	12.3031 mL	24.6063 mL
5 mM	0.4921 mL	2.4606 mL	4.9213 mL
10 mM	0.2461 mL	1.2303 mL	2.4606 mL
50 mM	0.0492 mL	0.2461 mL	0.4921 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481