


## Pefloxacin Mesylate Dihydrate

### Chemical Properties

CAS No. :	149676-40-4
Formula:	C <sub>17</sub> H <sub>20</sub> FN <sub>3</sub> O <sub>3</sub> ·CH <sub>4</sub> O <sub>3</sub> ·2H <sub>2</sub> O
Molecular Weight:	465.49
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



### Biological Description

Description	Pefloxacin Mesylate Dihydrate (1589 RB) , an antibacterial agent, restrains bacterial DNA replication by inhibiting DNA gyrase (topoisomerase).
Targets(IC50)	Antibacterial,Antibiotic,Parasite,DNA/RNA Synthesis,Topoisomerase
In vivo	Pefloxacin exerts potent bactericidal activity through several pathways and exhibits biphasic action: its killing effect decreases with increased concentration following initial rapid bactericidal activity. In bacterial cells cultured with Pefloxacin, the SOS response (non-replicative DNA repair) is induced, which inhibits replication, disrupts spindle formation, and blocks cell division, ultimately proving detrimental to bacteria by interfering with their morphological and biochemical properties. Similar to other quinolone antibiotics, Pefloxacin primarily targets bacterial DNA gyrase (topoisomerase II), a critical bacterial enzyme. It effectively inhibits members of the Enterobacteriaceae family (E. coli, Citrobacter, Klebsiella, Proteus, Morganella) with minimum inhibitory concentrations (MICs) ranging from 0.03 to 8 mg/L; for Shigella, Salmonella, and Yersinia, the MIC50/MIC90 values are 0.06/0.06, 0.12/1.0, and 0.12/0.25 mg/L respectively. Pefloxacin has weaker activity against Legionella pneumophila (MIC50/MIC90: 1.0/1.0 mg/L), very weak activity against Vibrio cholerae and Eikenella corrodens, and Listeria monocytogenes, Helicobacter pylori, and Nocardia asteroides show resistance with MIC50 ≥8 mg/L. Mycobacterial sensitivity to Pefloxacin varies from moderate to none, and most anaerobes are resistant. Rickettsiae and Chlamydiae are sensitive at treatable concentrations, while Mycoplasma and Ureaplasma show poor sensitivity (MIC50/MIC90: 2/8 mg/L). Pefloxacin's MICs for Gram-negative aerobes such as Alcaligenes, Pseudomonas aeruginosa, Acinetobacter, and Klebsiella range from 1 to 4 mg/L. For Gram-positive cocci, including coagulase-negative staphylococci and Staphylococcus aureus, even those resistant to other antibiotics, Pefloxacin maintains some sensitivity with MICs spanning 0.125 to 0.5 mg/L.

### Solubility Information

Solubility	DMSO: 1 mg/mL (2.15 mM),Sonication is recommended. Ethanol: < 1 mg/mL (insoluble or slightly soluble), H2O: 95 mg/mL (204.09 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

### Preparing Stock Solutions

---

	<b>1mg</b>	<b>5mg</b>	<b>10mg</b>
1 mM	2.1483 mL	10.7414 mL	21.4827 mL
5 mM	0.4297 mL	2.1483 mL	4.2965 mL
10 mM	0.2148 mL	1.0741 mL	2.1483 mL
50 mM	0.043 mL	0.2148 mL	0.4297 mL

---

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

#### Reference

Bergogne-Bérézin E, et al. Int J Antimicrob Agents. 1991;1(1):29-46.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481