

## Methyl-PEG4-acyl chloride

## Chemical Properties

CAS No. : 62124-69-0


Formula: C10H19ClO6

Molecular Weight: 270.71

Storage: Keep away from direct sunlight

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Methyl-PEG4-acyl chloride is a PEG-based linker for PROTACs that connects two essential ligands, facilitating the formation of PROTAC molecules and enabling selective protein degradation via the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands connected by a linker: one ligand targets an E3 ubiquitin ligase, while the other targets the protein of interest. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.694 mL	18.4699 mL	36.9399 mL
5 mM	0.7388 mL	3.694 mL	7.388 mL
10 mM	0.3694 mL	1.847 mL	3.694 mL
50 mM	0.0739 mL	0.3694 mL	0.7388 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481