

CRA-2059 TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	CRA-2059 is a highly specific and selective tryptase inhibitor with a K_i of 620 pM for recombinant human tryptase- β (rHT β) [1][2].
Targets(IC50)	Others,Serine Protease
In vitro	Tryptase, a trypsin-like serine protease, is a major protein in human mast cell secretory granules. CRA-2059 holds potential for inflammatory bowel disease research[1].

Reference

Tremaine WJ, et al. Treatment of mildly to moderately active ulcerative colitis with a tryptase inhibitor (APC 2059): an open-label pilot study. Aliment Pharmacol Ther. 2002;16(3):407-413.

Selwood T, et al. Potent bivalent inhibition of human tryptase-beta by a synthetic inhibitor. Biol Chem. 2003;384(12):1605-1611.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481