

## Monohydroxy Netupitant D6

## Chemical Properties


CAS No. :

Formula: C<sub>30</sub>H<sub>26</sub>D<sub>6</sub>F<sub>6</sub>N<sub>4</sub>O<sub>2</sub>

Molecular Weight: 600.63

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

Monohydroxy Netupitant-D6 is a metabolite of Netupitant (T12211).

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6649 mL	8.3246 mL	16.6492 mL
5 mM	0.333 mL	1.6649 mL	3.3298 mL
10 mM	0.1665 mL	0.8325 mL	1.6649 mL
50 mM	0.0333 mL	0.1665 mL	0.333 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481