


Corynecin I

Chemical Properties

CAS No. :	4423-58-9
Formula:	C ₁₁ H ₁₄ N ₂ O ₅
Molecular Weight:	254.242
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Corynecin I is a chloramphenicol-like antibiotic originally isolated from Corynebacterium. It is active against Gram-positive and Gram-negative bacteria with MIC values of 5.2-83 µg/ml for <i>S. sonnei</i> , <i>P. vulgaris</i> , <i>K. pneumoniae</i> , and <i>S. aureus</i> among others.
Targets(IC50)	Others,Antibiotic

Solubility Information

Solubility	DMSO: Soluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9333 mL	19.6665 mL	39.3329 mL
5 mM	0.7867 mL	3.9333 mL	7.8666 mL
10 mM	0.3933 mL	1.9666 mL	3.9333 mL
50 mM	0.0787 mL	0.3933 mL	0.7867 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481