

## Argiprestocin

## Chemical Properties

CAS No. : 9034-50-8

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Argiprestocin (Arginine vasotocin) , also called antidiuretic hormone that plays a key role in maintaining osmolality. It is a neurohypophysial hormone found in most mammals. Drugs used to cause constriction of the blood vessels.
Targets(IC50)	Others,Vasopressin Receptor

## Solubility Information

Solubility	H2O: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Reference

Tsuji T, et al. J Physiol. 2017 Apr 12.

Roussel R, et al. Curr Opin Nephrol Hypertens. 2017 Apr 11.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481