

Compound T128194

Chemical Properties


CAS No. :

Formula: C52H82O21

Molecular Weight: 1043.21

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2S,3R,4S,5S,6R)-6-(((2R,3R,4R,5S,6R)-3,4-dihydroxy-6-(hydroxymethyl)-5-(((2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy)oxan-2-yl]oxy)methyl)-3,4,5-trihydroxyoxan-2-yl (4aS,6aS,6bR,10S,12aR,14bS)-6a,6b,9,9,12a-pentamethyl-2-methylidene-10-(((2S,3R,4S,5S)-3,4,5-trihydroxyoxan-2-yl]oxy)-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,12b,13,14b-icosahydricene-4a-carboxylate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9586 mL	4.7929 mL	9.5858 mL
5 mM	0.1917 mL	0.9586 mL	1.9172 mL
10 mM	0.0959 mL	0.4793 mL	0.9586 mL
50 mM	0.0192 mL	0.0959 mL	0.1917 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481