


Filicane-3 $\beta$ ,4 $\alpha$ ,25-triol

## Chemical Properties

CAS No. : 2361548-00-5  
 Formula: C<sub>30</sub>H<sub>52</sub>O<sub>3</sub>  
 Molecular Weight: 460.743  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description Filicane-3 $\beta$ ,4 $\alpha$ ,25-triol is a natural product for research related to life sciences. The catalog number is TN5965 and the CAS number is 2361548-00-5.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1704 mL	10.8521 mL	21.7042 mL
5 mM	0.4341 mL	2.1704 mL	4.3408 mL
10 mM	0.217 mL	1.0852 mL	2.1704 mL
50 mM	0.0434 mL	0.217 mL	0.4341 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481