


MK190

## Chemical Properties

CAS No. : 4204-99-3  
Formula: C<sub>10</sub>H<sub>8</sub>FN<sub>3</sub>O<sub>2</sub>  
Molecular Weight: 221.19  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	MK190 is a potential bioactive agent. No further details yet.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.521 mL	22.605 mL	45.210 mL
5 mM	0.9042 mL	4.521 mL	9.042 mL
10 mM	0.4521 mL	2.2605 mL	4.521 mL
50 mM	0.0904 mL	0.4521 mL	0.9042 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481