


GLP-1R agonist 10

Chemical Properties

CAS No. :	2762515-25-1
Formula:	C ₂₉ H ₂₉ ClFN ₅ O ₄
Molecular Weight:	566.02
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	GLP-1R agonist 10 is a GLP-1 agonist (EC ₅₀ : 0.051 nM).
Targets(IC ₅₀)	Others, Glucagon Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7667 mL	8.8336 mL	17.6672 mL
5 mM	0.3533 mL	1.7667 mL	3.5334 mL
10 mM	0.1767 mL	0.8834 mL	1.7667 mL
50 mM	0.0353 mL	0.1767 mL	0.3533 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481