

Anti-PIKFyve

Chemical Properties

CAS No. :

Formula: C₁₀H₁₁F₂N₂Na₄O₁₃P₃S

Molecular Weight: 622.14

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	antibody against human PIKfyve
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6074 mL	8.0368 mL	16.0736 mL
5 mM	0.3215 mL	1.6074 mL	3.2147 mL
10 mM	0.1607 mL	0.8037 mL	1.6074 mL
50 mM	0.0321 mL	0.1607 mL	0.3215 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481