


## BNP (1-32), human

## Chemical Properties

CAS No. :	114471-18-0
Formula:	C143H244N50O42S4
Molecular Weight:	3464.04
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	BNP (1-32), human is a brain natriuretic peptide
Targets(IC50)	RAAS,Others

## Solubility Information

Solubility	DMSO: 206.6 mg/mL (59.64 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.2887 mL	1.4434 mL	2.8868 mL
5 mM	0.0577 mL	0.2887 mL	0.5774 mL
10 mM	0.0289 mL	0.1443 mL	0.2887 mL
50 mM	0.0058 mL	0.0289 mL	0.0577 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Tsuruda T , Boerrigter G , Huntley B K , et al. Brain Natriuretic Peptide Is Produced in Cardiac Fibroblasts and Induces Matrix Metalloproteinases[J]. Circulation Research, 2002, 91(12):1127.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481