

## 2-acetamido-3-(3-hydroxy-4-methoxy-5-methylphenyl)propanoic acid

## Chemical Properties


CAS No. :

Formula: C<sub>13</sub>H<sub>17</sub>NO<sub>5</sub>

Molecular Weight: 267.28

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-acetamido-3-(3-hydroxy-4-methoxy-5-methylphenyl)propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7414 mL	18.707 mL	37.4139 mL
5 mM	0.7483 mL	3.7414 mL	7.4828 mL
10 mM	0.3741 mL	1.8707 mL	3.7414 mL
50 mM	0.0748 mL	0.3741 mL	0.7483 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481